



Novel Visual Feedback Device Provides Enhanced Patient Coaching Experience for Patients Undergoing Left Breast DIBH

Presenting Author: David D. Campos

Session Title: Treatment Implications of Motion Management

Date: 8/2/2018


Time: 7:30-9:30AM



Sample breath-holds from two patients on trial. Breathe Well agrees temporally with AlignRT breathing traces.



- couch-mounted
- vendor-independent
- adjustable settings (A 1,2,3)
- optical and depth information (B)
- real-time feedback via patient screen



Chest-wall position data from two patients on trial shows **significant improvement in reproducibility** when using visual feedback for coaching the patients